

# ScleroSafe™



The ScleroSafe platform is an advanced non-thermal, non-tumescent technology solution that eliminates the need for surgery or costly capital equipment, while minimizing the pain and recovery time associated with conventional treatment.

A minimally invasive, tumescent-free approach to simply and quickly eliminate varicose veins and deliver immediate, long-lasting results by ablating superficial veins with a diameter of 2 to 5mm.

ScleroSafe can be used with the 5cc\10cc Dual Syringe, a specially designed 5F 150mm long or 5F 350mm long double lumen catheter and a dedicated guide wire.



ScleroSafe inverse action dual 10 ml / 5 ml syringe



ScleroSafe double lumen catheter  
(5F/ 150mm or 5F/ 350mm)

Order Code	Description	Master Pack
VAR0150	ScleroSafe - 150mm	10 units
VAR0350	ScleroSafe - 350mm	10 units